


Committed to excellence



RUTRONIK
ELECTRONICS WORLDWIDE



FFP2 FACE MASKS & SARS-COV-2 ANTIGEN RAPID TEST

corona@rutronik.com

KRÖNER
MEDIZINTECHNIK
Technik und Leidenschaft



SARS-CoV-2 Antigen Rapid Test

The Rapid Antigen Test for Novel Coronavirus (SARS-Cov-2) is an in vitro diagnostic test for the qualitative detection of novel coronavirus antigens in nasopharyngeal swabs and oropharyngeal swabs using the rapid immunochromatographic method. Identification is based on the monoclonal antibodies specific for the novel coronavirus antigen. This product is intended for professional use in the laboratory and at the point of care.

Features

- High clinical sensitivity (92.00%) and specificity (99.26%)
- Rapid and reliable test results in 15 minutes
- Can be performed with nasal samples
- All test components are included
- Tested by Paul-Ehrlich-Institute
- Officially listed by the federal institute for pharmaceuticals and medical products




available at Rutronik24.com




PU 25
Test carrier, package insert, sterilised swab, nozzle with filter and sample collection buffer


Operation steps




Step1: Collect the swab



Step2: Add 4 drops of buffer to hole A




Step3: Insert the sample through hole B to hole A, rotate the swab clockwise first and counterclockwise, twice each.



Step4: Peel off the adhesive, and fold the left side over, wait for 15 minutes.


Result Interpretation

Negative




Top line only

Positive



Both bottom line and top line

Invalid



Bottom line only, or no line

The current test card is based on the specific antibody-antigen reaction and immunoanalysis technology. The test card contains colloidal gold labeled SARS-CoV-2 N protein monoclonal antibody which is pre-coated on the combination pad, matched SARS-CoV-2 N protein monoclonal antibody immobilized on the Test area and corresponding antibody in the quality control area. During testing, the N protein in the sample combines with the colloidal gold labeled SARS-CoV-2 N protein monoclonal antibody which is pre-coated on the combination pad. The conjugates migrate upward under capillary effect, and subsequently captured by the N protein monoclonal antibody immobilized in the Test area. The higher the contents of N protein in the sample, the more the conjugates captures and the darker the color in the test area is. If there is no virus in the sample or the virus content is lower than the detection limit, then there is no color demonstrated in the test area. Regardless of the presence or absence of the virus in the sample, a purple stripe will appear in the quality control area. The purple stripe in the quality control area is a criterion for the judgment of whether or not there is enough sample and whether or not the chromatography procedure is normal.

Contact: **Rutronik Elektronische Bauelemente GmbH** | Tilo Rollwa Director Digital Marketing - Tel.: +49 7231 801 15 40 | tilo.rollwa@rutronik.com
 Stefanie Piller - Tel.: +49 7231 801 15 42 | stefanie.piller@rutronik.com | Nikolai Schnarz - Tel.: +49 7231 801 14 27 | nikolai.schnarz@rutronik.com



Document No.: CE-DOC-CG27
Rev.: 1/0

Declaration of Conformity

Manufacture Address: Beijing Lepu Medical Technology Co., Ltd.
Building 7-1 No.37 Chaoqian Road, Changping District,
Beijing, 102200, P.R. China

European Representative: Lepu Medical (Europe) Cooperatief U.A.
Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, The
Netherlands

Product information: SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold
Immunochromatography)
Model:
1 test/kit; 5 tests/kit; 10 tests/kit; 25 tests/kit; 50 tests/kit

Classification: Others (not in List A and List B)

Conformity Assessment Route: Section 2 to 5 in annex III of IVDD 98/79/EC
We herewith declare that the above mentioned products
meet the provisions of the following EC Council Directives
and Standards.
All supporting documentations are retained under the
premise of the manufacturer.

General Applicable Directive: DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT
AND OF THE COUNCIL of 27 October 1998 on *in vitro*
diagnostic medical devices

Standards Applied: All applicable harmonized standards (published in the
official journal of the European Communities on 25th March
2020).
The applicable standards are listed in Annex 1.

Place, date of issue Beijing, P.R. China, 3th, Sept., 2020

**Signature of Management
Representative**

Zhao Anqun

Beijing Lepu Medical Technology Co., Ltd.
Building 7-1 No.37 Chaoqian Road, Changping District, Beijing, 102200, P.R. China



Annex 1

EN ISO 13485:2016 Medical devices – quality management systems - requirements
for regulatory purposes
EN ISO 14971:2019 Medical devices – application of risk management to medical
devices
EN ISO 15223-1:2016 Medical devices – Symbols to be used with medical device
labels, labelling and information to be supplied – Part 1: General requirements
EN 1041:2008+A1:2013 Information supplied by the manufacturer of medical devices
EN ISO 18113-1:2011 In vitro diagnostic medical devices – information supplied by
the manufacturer (labelling) – Part 1: terms, definitions and general requirements
EN ISO 18113-2:2011 In vitro diagnostic medical devices – information supplied by
the manufacturer (labelling) – Part 2: in vitro diagnostic reagents for professional use
EN ISO 23640:2015 In vitro diagnostic medical devices – evaluation of stability of in
vitro diagnostic reagents
EN 13612:2002/AC: 2002 Performance evaluation of in vitro diagnostic medical
devices
IEC 62366-1:2015 Application of usability engineering to medical devices

FFP2
EN149:2001+A1:2009 **SANOCARE+**

KRÖNER
MEDIZINTECHNIK
Technik und Leidenschaft

Particle Filtering Face Mask

The **SANOCARE+** FFP2 mask offers a comfortable fit thanks to the pleasantly soft inner fleece and extra-wide, tight-fitting ear straps. It meets all standards and certifications. The 5 layers offer optimal protection against fine dust and aerosols. The flexible nasal clip enables a good and firm fit.




Brand	SANOCARE+
Model #	FH-1
PU	20 pcs (single packaging)
Certification #	CE 0598
Tested according to norm	EN 149:2001+A1:2009
Filter	FFP2
Layers	5
Nasal clip	Yes
Material	54% Fleece 24,5% Meltblown- Fleece 21,5% Hot-Air Cotton
Reusable	No


available at [Rutronik24.com](https://www.rutronik24.com)



CE 0598



ANLEGEINSTRUKTION



1 Maske mit Hilfe der Bänder auf falten. Außenseite mit Aufdruck muss nach vorne zeigen, der Nasenbügel befindet sich oben.



2 Maske über Mund und Nase positionieren, beide Gummibänder hinter den Kopf führen und Schlaufen hinter den Ohren fixieren.



3 Den Nasenbügel mithilfe beider Hände durch leichten Druck an die Nasenkontur angleichen für einen bestmöglichen Dichtsitz.



4 Dichtsitz kontrollieren: Beide Hände um die Maske legen, ohne dass diese sich dabei bewegt, und kräftig ausatmen. Wenn Luft um Ihre Nase entweicht, drücken Sie den Bügel fester um Ihren Nasenrücken.

WARNUNGEN & HINWEISE

- Die Warnungen & Hinweise müssen vor Verwendung der Maske gelesen und verstanden werden. Falsche Anwendung oder Nichtbeachtung kann zu Krankheit, Verletzungen und Gesundheitsschäden führen. Maske vor Gebrauch prüfen.
- Stellen Sie sicher, dass die Klassifizierung der Maske für den vorgesehenen Einsatzzweck geeignet ist. Der Anwender trägt die Verantwortung dafür, dass die Atemschutzmaske die notwendige Schutzstufe für die Art und Konzentration der Verschmutzung in dem Bereich, in dem sie angewendet werden soll, erreicht.
- Die Maske muss dicht sitzen. Führen Sie vor Verwendung unbedingt eine Dichtsitzkontrolle durch.
- Es ist unwahrscheinlich, dass ein ausreichender Dichtsitz erreicht wird, falls sich Gesichtshaare wie Bärte oder Koteletten unter der Dichtlinie am Gesicht befinden.
- Mindestens 18 Vol.% Sauerstoffgehalt in der Atemluft sind notwendig, um die Maske sicher benutzen zu können.

- Diese Maske mit der Kennzeichnung NR ist zum einmaligen Gebrauch bestimmt und ist nach einer Schicht bzw. maximal 8h Gebrauchsdauer zu vernichten.
- Bei Beschädigungen der Maske oder bei ansteigendem Atemwiderstand ist die Maske umgehend abzulegen und zu wechseln.
- Die Maske sollte in Innenräumen bei einer Temperatur von -25 bis 40°C sowie einer Luftfeuchtigkeit von maximal 80% gelagert werden. Direkte Sonneneinstrahlung und Feuchtigkeit sind zu vermeiden.
- Nehmen Sie keine Veränderungen oder Reparaturen an der Maske vor.
- Die Maske bietet Schutz vor festen und flüssigen Aerosolen im Rahmen der Schutzklasse FFP2, kann jedoch die Möglichkeit einer Ansteckung oder Erkrankung nicht gänzlich eliminieren.

EINSATZBEREICH

Diese Atemschutzmaske bietet Schutz vor festen und flüssigen Aerosolen bis zum 10-fachen des Arbeitsplatzgrenzwertes sowie zusätzlich gegen biologische Arbeitsstoffe der Gruppe 2.

MATERIAL & HALTBARKEIT

54% Vlies, 24,5% Schmelzgeblasener Vliesstoff, 21,5% Heißluftbaumwolle

- Haltbarkeit: 3 Jahre
- Herstellungsdatum: siehe Beipackzettel
- Chargennummer: siehe Beipackzettel

Hersteller: Name, Anschrift
Shenzhen Baishi Health Medicine Co.,Ltd, 601, 10TH Plant, Longxing Street, Dakang Community, Yuan shan road, Long gang District, Shenzhen city, P.R.China

Beauftragte notifizierte Stelle der Konformitätsbewertung: Name, Anschrift, Kennnummer
SGS FIMKO OY, Takomitie 8, FI-00380 Helsinki, Finland, Tel: +358 9 696 361, www.sgs.com / CE 0598

Link zur EU Konformitätserklärung zu finden unter:
www.hms24.eu/sanocare-plus/

Name und Fundstelle der Verordnung:
Verordnung (EU) 2016/425, Amtsblatt der Europäischen Union, Nr. L 81/51, 2016

Name und Fundstelle der Norm:
EN 149:2001+A1:2009 www.beuth.de





FFP2 NR
geprüft nach EN 149:2001+A1:2009

SANOCARE+
Partikelfiltrierende Atemschutzmaske

- Herstellerinformationen sind unbedingt zu beachten
- Zulässiger Temperaturbereich bei Lagerung: -25°C bis +40°C
- Zulässige relative Luftfeuchte bei Lagerung: max. 80%
- Zur einmaligen Verwendung



MADE IN CHINA



SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 2 of 10
 Test Result

Personal Protective Equipment - Respiratory Protective Devices- Filtering Half Masks to Protect against Particles- Requirements, Testing, Marking
 EN 149:2001+A1:2009

Clause 7.4 Packaging
 (EN 149:2001+A1:2009 Clause 8.2)

Test Requirement	Results	Comment
Particle filtering half masks shall be offered for sale packaged in such a way that they are protected against mechanical damage and contamination before use.	Comply	Pass

Clause 7.5 Material
 (EN 149:2001+A1:2009, Clause 8.2 & 8.3.1 & 8.3.2)

Test Requirement	Results	Comment
Materials used shall be suitable to withstand handling and wear over the period for which the particle filtering half mask is designed to be used.	Comply	Pass
After undergoing the conditioning described in 8.3.1 none of the particle filtering half masks shall have suffered mechanical failure of the facepiece or straps.	Comply	
When conditioned in accordance with 8.3.1 and 8.3.2 the particle filtering half mask shall not collapse.	Comply	
Any material from the filter media released by the air flow through the filter shall not constitute a hazard or nuisance for the wearer.	Comply	

Clause 7.6 Cleaning and Disinfecting
 (EN 149:2001+A1:2009, Clause 8.4 & 8.5 & 8.11)

Test Requirement	Results	Comment
If the particle filtering half mask is designed to be re-usable, the materials used shall withstand the cleaning and disinfecting agents and procedures to be specified by the manufacturer.	Not applicable (Not designed to be re-usable)	N.A.
With reference to 7.9.2, after cleaning and disinfecting the re-usable particle filtering half mask shall satisfy the penetration requirement of the relevant class.		

Clause 7.7 Practical Performance
 (EN 149:2001+A1:2009, Clause 8.4)

Test Requirement	Results	Comment
The particle filtering half mask shall undergo practical performance tests under realistic conditions. These general tests serve the purpose of checking the equipment for imperfections that cannot be determined by the tests described elsewhere in this standard.	No imperfections	Pass

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 3 of 10
 Test Result

Clause 7.8 Finish of Parts
 (EN 149:2001+A1:2009, Clause 8.2)

Test Requirement	Results	Comment
Parts of the device likely to come into contact with the wearer shall have no sharp edges or burrs.	Comply	Pass

Clause 7.9.1 Total Inward Leakage
 (EN 149:2001+A1:2009, Clause 8.5)

Test Requirement	Results	Comment
The total inward leakage consists of three components: face seal leakage, exhalation valve leakage (if exhalation valve fitted) and filter penetration. For particle filtering half masks fitted in accordance with the manufacturer's information, at least 46 out of the 50 individual exercise results (i.e. 10 subjects x 5 exercises) for total inward leakage shall be not greater than: 25% for FFP1, 11% for FFP2, 5% for FFP3	Detail refer to Appendix 1	Meet FFP1, Meet FFP2

Appendix 1: Summarization of Test Data

Inward Leakage Test Data

Subject	Sample No.	Condition	Walk (%)	Head Side/side (%)	Head up/down (%)	Talk (%)	Walk (%)	Mean (%)
Zhou	1	A.R.	7.00	6.77	5.84	5.10	6.91	6.32
Luo	2	A.R.	5.82	6.49	6.46	7.53	7.33	7.13
Lu	3	A.R.	5.36	5.60	7.45	7.08	6.40	6.38
Wang	4	A.R.	5.22	4.65	5.79	6.08	4.91	5.33
Bao	5	A.R.	7.72	7.88	7.32	6.49	6.77	7.24
Ding	6	T.C.	5.81	5.34	6.36	4.77	4.50	5.44
Li	7	T.C.	7.53	6.56	7.52	6.48	8.70	7.76
Chen	8	T.C.	6.36	4.57	5.84	5.41	4.88	5.41
Song	9	T.C.	7.29	5.85	5.59	7.04	6.46	6.45
Ye	10	T.C.	6.85	7.18	7.53	6.67	8.94	7.83

Facial Dimension (mm)

Subject	Face length	Face Width	Face Depth	Mouth Width
Chen	125	150	120	58
Lu	115	132	107	48
Zhou	115	136	106	52
Li	126	130	107	46
Luo	125	136	100	43
Zheng	128	140	112	55
Wang	120	147	103	48
Song	120	140	100	50
Bao	130	134	104	50

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 6 of 10
 Test Result

Clause 7.13 Head Harness
 (EN 149:2001+A1:2009, Clause 8.4 & 8.5)

Test Requirement	Results	Comment
The head harness shall be designed so that the particle filtering half mask can be donned and removed easily.	Comply	Pass
The head harness shall be adjustable or self-adjusting and shall be sufficiently robust to hold the particle filtering half mask firmly in position and be capable of maintaining total inward leakage requirements for the device.	Comply	

Clause 7.14 Field of Vision
 (EN 149:2001+A1:2009, Clause 8.4)

Test Requirement	Results	Comment
The field of vision is acceptable if determined so in practical performance tests.	Comply	Pass

Clause 7.15 Exhalation Valve(s)
 (EN 149:2001+A1:2009, Clause 8.2 & 8.9.1 & 8.3.4 & 8.8)

Test Requirement	Results	Comment
(a) A particle filtering half mask may have one or more exhalation valve(s), which shall function correctly in all orientations.	Not applicable due to No exhalation valve	N.A.
(b) If an exhalation valve is provided it shall be protected against or be resistant to dirt and mechanical damage and may be shrouded or may include any other device that may be necessary for the particle filtering half mask to comply with 7.9.	Not applicable due to No exhalation valve	
(c) Exhalation valve(s), if fitted, shall continue to operate correctly after a continuous exhalation flow of 300 l/min over a period of 30 s.	Not applicable due to No exhalation valve	
(d) When the exhalation valve housing is attached to the face/bank, it shall withstand axially a tensile force of 10N applied for 10 s.	Not applicable due to No exhalation valve	

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 7 of 10
 Test Result

Clause 7.16 Breathing Resistance
 (EN 149:2001+A1:2009, Clause 8.9)

Test Requirement	Results	Comment
The penetration of the filter of the particle filtering half mask shall meet the requirements of the following table.	Detail refer to Appendix 5	Meet FFP1, Meet FFP2, Meet FFP3

Appendix 5: Summarization of Test Data

Breathing resistance (mbar)

Classification	Maximum permitted resistance (mbar)	1						2						3								
		Inhalation	30 l/min		95 l/min		Exhalation	30 l/min		95 l/min		Inhalation	30 l/min		95 l/min		Exhalation	30 l/min		95 l/min		
FFP1	0.6	0.6	2.1	2.1	3.0	3.0	0.6	2.1	2.1	3.0	0.6	2.1	2.1	3.0	0.6	2.1	2.1	3.0	0.6	2.1	2.1	3.0
FFP2	0.7	0.7	2.4	2.4	3.0	3.0	0.7	2.4	2.4	3.0	0.7	2.4	2.4	3.0	0.7	2.4	2.4	3.0	0.7	2.4	2.4	3.0
FFP3	1.0	1.0	3.0	3.0	3.0	3.0	1.0	3.0	3.0	3.0	1.0	3.0	3.0	3.0	1.0	3.0	3.0	3.0	1.0	3.0	3.0	3.0

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 4 of 10
 Test Result

Classification	Sample No.	Condition	Walk (%)	Head Side/side (%)	Head up/down (%)	Talk (%)	Walk (%)	Mean (%)
Ding	134	150	110	52				
Liu	120	135	117	50				
Ye	126	137	105	52				

Clause 7.9.2 Penetration of Filter Material
 (EN 149:2001+A1:2009, Clause 8.11 & EN 13274-7:2019)

Test Requirement	Results	Comment
The penetration of the filter of the particle filtering half mask shall meet the requirements of the following table.	Detail refer to Appendix 2	Meet FFP1, Meet FFP2

Appendix 2: Summarization of Test Data

Penetration of filter material

Aerosol	Condition	Sample No.	Penetration (%)
Sodium chloride test	As received	1	0.269
		2	0.274
		3	0.263
	Simulated wearing treatment	4	0.265
		5	0.274
		6	0.283
Mechanical strength + Temperature conditioned	7	0.347	
	8	0.425	
	9	0.365	
Paraffin oil test	As received	10	1.609
		11	1.556
		12	1.527
	Simulated wearing treatment	13	1.561
		14	1.630
		15	1.704
Mechanical strength + Temperature conditioned	16	4.681	
	17	5.237	
	18	5.014	

Flow conditioning : Single filter: 95.0 L/min

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 5 of 10
 Test Result

Clause 7.10 Compatibility with Skin
 (EN 149:2001+A1:2009, Clause 8.4 & 8.5)

Test Requirement	Results	Comment
Materials that may come into contact with the wearer's skin shall not be known to be likely to cause irritation or any other adverse effect to health.	No irritation or any other adverse effect to health	Pass

Clause 7.11 Flammability
 (EN 149:2001+A1:2009, Clause 8.6)

Test Requirement	Results	Comment
The material used shall not present a danger for the wearer and shall not be of highly flammable nature	Detail refer to Appendix 3	Pass
When tested, the particle filtering half mask shall not burn or not to continue to burn for more than 5 s after removal from the flame.		

Appendix 3: Summarization of Test Data

Flammability

Condition	Sample No.	Result
As received	1	NIL
	2	NIL
	3	NIL
Temperature conditioned	4	NIL

Clause 7.12 Carbon Dioxide Content of The Inhalation Air
 (EN 149:2001+A1:2009, Clause 8.7)

Test Requirement	Results	Comment
The carbon dioxide content of the inhalation air (dead space) shall not exceed an average of 1.0 % (by volume)	Detail refer to Appendix 4	Pass

Appendix 4: Summarization of Test Data

Carbon Dioxide Content of The Inhalation Air

Condition	Sample No.	Result (%)
As received	1	0.4637
	2	0.4627
	3	0.4630

Mean value: 0.46

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 8 of 10
 Test Result

Clause 7.17 Clogging
 (EN 149:2001+A1:2009, Clause 8.9 & 8.10)

Test Requirement	Results	Comment
Clause 7.17.2 Breathing resistance Valved particle filtering half masks: After clogging the inhalation resistances shall not exceed: FFP1: 4 mbar, FFP2: 5 mbar, FFP3: 7 mbar at 95 l/min continuous flow The exhalation resistance shall not exceed 3 mbar at 160 L/min continuous flow.	Optional for single shift device only	N.A.
Valveless particle filtering half masks: After clogging the inhalation and exhalation resistances shall not exceed: FFP1: 3 mbar, FFP2: 4 mbar, FFP3: 5 mbar at 95 l/min continuous flow		
Clause 7.17.3 Penetration of filter material All types (valved and valveless) of particle filtering half masks claimed to meet the clogging requirement shall also meet the requirements.	Optional for single shift device only	N.A.

Classification	n	Maximum penetration of test aerosol	%
FFP1	20	Sodium chloride test 95 l/min	20
		Paraffin oil test 95 l/min	6
		% max.	1
FFP2	6	Sodium chloride test 95 l/min	20
		Paraffin oil test 95 l/min	6
		% max.	1
FFP3	1	Sodium chloride test 95 l/min	20
		Paraffin oil test 95 l/min	6
		% max.	1

Clause 7.18 Detachable Parts
 (EN 149:2001+A1:2009, Clause 8.2)

Test Requirement	Results	Comment
All detachable parts (if fitted) shall be readily connected and secured, where possible by hand	No detachable parts	N.A.

Test	Uncertainty
Total inward leakage	3.4%
Penetration of filter material	4.8%
Carbon dioxide content of the inhalation air	3.9%
Breathing resistance (30L/min)	5.9%
Breathing resistance (95L/min)	4.9%
Breathing resistance (160L/min)	4.3%

Member of the SGS Group (SGS SA)

SGS
 Test Report SL52025256705901TX Date: June 13, 2020 Page 9 of 10
 Test Result

Sample Photo

Member of the SGS Group (SGS SA)

Price List (Net unit price in €)

PU 25 pcs.



**DISMED1002
Corona Antigen Rapid Test**

25 Pcs.*	6,95/Pcs.
100 Pcs.	4,95/Pcs.
500 Pcs.	4,45/Pcs.
1.000 Pcs.	3,95/Pcs.
10.000 Pcs.	3,25/Pcs.

*MOQ

PU 20 pcs.



DISMED1001 FFP2-Mask

100 Pcs.*	0,65/Pcs.
1.000 Pcs.	0,59/Pcs.
10.000 Pcs.	0,55/Pcs.

*MOQ